



## **NO BIOHAZARDS**

**DO NOT** use this autoclave for the decontamination of biological and potentially infectious waste. Sterilization of glassware and liquids may be done in this particular autoclave.

- ◆ Never put sealed containers, solvents, volatile or corrosive chemicals, or radioactive materials in the autoclave.
- ◆ If the cycle did not meet the appropriate operating conditions or aborted, repeat the autoclave cycle. If the cycle fails again, report the malfunction to the responsible person listed below.
- ◆ Before unloading the Autoclave, wait for the pressure gauge to drop to zero with zero time remaining before opening the door.
- ◆ Using proper heat resistant PPE, carefully remove autoclaved materials after the cycle is completed and the appropriate waiting time for liquids has been observed.

### **Responsible Person for this Autoclave**

Name: \_\_\_\_\_

Phone: \_\_\_\_\_



<http://ehs.uky.edu/biosafety/>  
[biosafety@uky.edu](mailto:biosafety@uky.edu)